# MRI stimulus protocol

April 4, 2023

## 1 MonBinMel MRI design

This project will examine how and where ipRGC photoreceptor signals from each eye are integrated in the brain. The proposed design combines square wave monocular and binocular melanopsin-directed contrast transitions with brief sinusoidal LMS probes of varying intensity. This will allow us to map BOLD activation for monocular and binocular melanopsin stimulation in a basic block design whilst also measuring behavioral and BOLD responses to LMS probes under different levels of melanopsin stimulation.

To date, only a handful of fMRI studies have examined the BOLD response to melanopsin stimulation. Spitschan et al., (2016) measured BOLD activation to melanopsin-directed flicker across a range of frequencies (0.5, 1.0, 2.0, 4.0, 8.0, 16.0, 32.0, 64.0 Hz) and did not find a cortical response that exceeded that of the control conditions. In a subsequent study, Spitschan et al. (2017) presented 3 s pulses of melanopsin-directed contrast (25%, 50%, 100%, 200%, 400%) and observed BOLD activation in primary visual cortex that scaled with the contrast. The subjects in this experiment also reported that high-melanopsin contrast stimuli elicited a distinct perceptual experience, summarily described as an "unpleasant, blurry, minimal brightening that quickly faded".

Both of the aforementioned studies involved monocular stimulation with low sample sizes ( $n \le 4$ ) and focused on cortical visual areas (e.g., V1-3, LOC, MT) in the analyses. However, a recent analysis by Spitschan and Cajochen (2019) of existing published data implies a substantial binocular facilitation effect for melanopsin stimulation with respect to melatonin suppression. This putative effect may present as increased BOLD activation in subthalamic ipRGC project sites such as the geniculate, suprachiasmatic and pretectal nuclei. The current design permits a LowMel < MonMel < BinMel analysis and is ideally suited to exploring the BOLD representation of this effect. Further, there is reason to believe that melanopsin contrast facilitates the processing of cone signals and influences brightness perception (e.g., Brown et al., 2012; Allen et al., 2014; Davis et al., 2015), so we might reasonably expect improved detection for the LMS probes and a greater BOLD signal in primary visual areas during the high-melanopsin conditions.

### Hypotheses

- $H_{0,1}$  The BOLD signal in ipRGC projection targets (LGN, SCN, PN) will not differ across levels of melanopic stimulation.
- $H_{1,1}$  The BOLD signal in ipRGC projection targets (LGN, SCN, PN) will scale with the levels of melanopic stimulation (LowMel < MonMel < BinMel).
- $H_{0,2}$  Behavioural detection of LMS probes and associated cortical BOLD signals will not differ across melanopsin conditions.
- $H_{1,2}$  Behavioural detection of LMS probes and associated cortical BOLD signals will be facilitated by higher levels of melanopic stimulation.

#### 1.1 Basic block design

Each block will last 360 s (TR = 3 s) and employ a basic 30 s off (LowMel) / 30 s on (HighMel) design. The HighMel stimulation is subdivided into monocular and binocular conditions. We therefore have the following conditions:

- x6 LowMel (low melanopsin spectrum in both eyes)
- x2 MonLMel (high melanopsin spectrum in the left eye)
- x2 MonRMel (high melanopsin spectrum in the right eye)
- x2 **BinMel** (high melanopsin spectrum in both eyes)

Each 30 s period will include four LMS probes at random non-overlapping time points. These probes serve as attention events but will also allow us to measure contrast sensitivity functions separately for monocular and binocular melanopsin conditions. The probes are 4 Hz sinusoidal LMS contrast modulations lasting 1 s, and will not occur within 3 s of a spectral transition. The following (michelson) contrast levels should be achievable:

- 10.0%
- 5.0%
- 2.5%
- 1.25%

The code below creates the basic block design described above.


```
[78]: import random
      import numpy as np
      import matplotlib.pyplot as plt
      plt.style.use('bmh')
      from pysilsub import waves
      # https://www.geeksforgeeks.org/python-non-overlapping-random-ranges/
      def non_overlapping_ranges(n, width, low, high):
          """Find non-overlapping ranges for LMS probes."""
          # Maybe add an additional spacer parameter
          ranges = set()
          for _ in range(n):
              candidate = random.randint(low, high - width)
              while any(candidate >= (idx-(width*2))
                        and candidate <= (idx + (width*2))</pre>
                        for idx in ranges):
                  candidate = random.randint(0, high - width)
              ranges.add(candidate)
          ranges = [(idx, idx + width) for idx in ranges]
          return ranges
      # Parameters
      Fs = 50 # Sampling frequency of STLABs
```

```
block_length = 360000 # ms
duration = 30  # Basic on-off duration
# Basic square wave
off = np.zeros(duration*Fs)
on = np.ones(duration*Fs)
# Full cycle for mel-on and mel-off conditions
meloff = np.hstack([off, off])
melon = np.hstack([off, on])
# Make conditions (separate stimulus profiles for each eye)
conditions = ['MonLMel', 'MonRMel', 'BinMel']
left_stims = []
right_stims = []
for i, c in enumerate(conditions * 2):
   if c == 'MonLMel':
       left = melon.copy()
       right = meloff.copy()
   if c == 'MonRMel':
       left = meloff.copy()
       right = melon.copy()
   if c == 'BinMel':
       left = melon.copy()
       right = melon.copy()
   # Save the stims
   left_stims.append(left)
   right_stims.append(right)
# Shuffle up the stims
shuffled = [(1, r, c) for (1, r, c) in zip(left_stims, right_stims, conditions_
→* 2)]
random.shuffle(shuffled)
left_stims = [s[0] for s in shuffled]
right_stims = [s[1] for s in shuffled]
cons = [s[2] for s in shuffled]
# Full block melanopsin profile for each eye
left_mel = np.hstack(left_stims)
right_mel = np.hstack(right_stims)
time = np.linspace(0, block_length, len(right_mel)) / 1000
# Full block LMS profile (same for each eye)
```

```
lms mod = waves.make_stimulus_waveform(frequency=4, sampling_frequency=Fs,__
 →duration=1)
lms = np.ones(len(left_mel)) / 2
lms conditions = [.1, .05, .025, .0125]
lms contrasts = []
lms idxs = []
for idx in range(0, len(right_mel), 1500):
    # Use the non-overlapping ranges function to get idxs for LMS mods
    ranges = non_overlapping_ranges(n=4, width=50, low=150, high=1350)
    new = [(start+idx, end+idx) for (start, end) in ranges]
    lms idxs += new
    random.shuffle(lms_conditions)
    lms_contrasts += lms_conditions
for i, c in zip(lms_idxs, lms_contrasts):
    lms[i[0]:i[1]] = lms[i[0]:i[1]] + (lms_mod*c*5)
# Pl.ot.
fig, axs = plt.subplots(2, 1, figsize=(12, 4))
axs[0].plot(time, right mel, lw=3, c='k')
axs[0].set ylabel('Right eye\nMelanopsin')
axs[1].plot(time, left_mel, lw=3, c='k')
axs[1].set_ylabel('Left eye\nMelanopsin')
axs[1].set_xlabel('Time (s)')
for ax in axs:
    ax.set_ylim((-.5, 1.5))
    ax.set_yticks((0,1), ['Low', 'High'])
# lms axis
ax1 = axs[0].twinx()
ax2 = axs[1].twinx()
ax1.plot(time, lms, lw=1, c='b')
ax2.plot(time, lms, lw=1, c='b')
ax1.set_ylabel('Right eye\nLMS', color='b')
ax2.set_ylabel('Left eye\nLMS', color='b')
for ax in [ax1, ax2]:
    ax.set_ylim((-.5, 1.5))
    ax.tick_params(axis='y', colors='b')
    ax.set_yticks((0,1), ['-10\%', '10\%'])
    ax.set_xticks(np.arange(min(time), max(time)+1, 60))
# Add conditions
for x, t in zip(range(45, 361, 60), cons):
    plt.text(x, 1.8, t, horizontalalignment='center', fontsize=12)
```

```
for x, t in zip(range(15, 361, 60), ['LowMel']*6):
    plt.text(x, 1.8, t, horizontalalignment='center', fontsize=12)

fig.suptitle('Block design - melanopsin with LMS contrast probes')
plt.tight_layout()
```



```
[79]: # Save the events for FSL
     events = (
         pd.concat(
                pd.DataFrame(zip(range(30, 360, 60), cons), columns=['Onset',__
      pd.DataFrame(zip(range(0, 360, 60), ['LowMel']*6), __

columns=['Onset', 'Event']),
                pd.DataFrame(sorted(zip([1[0]/50 for 1 in lms_idxs],__
      columns=['Onset', 'Event', 'Cm'])
            ]
         ).sort_values(by='Onset').reset_index(drop=True)
     pd.set option('display.max rows', 100)
     events.loc[events['Event']=='LMS', 'Duration'] = 1
     events.loc[events['Event']!='LMS', 'Duration'] = 30
     events
```

```
[79]:
           Onset
                     Event
                                 Cm Duration
      0
            0.00
                    LowMel
                                NaN
                                          30.0
      1
            7.18
                       LMS
                             0.0250
                                           1.0
      2
           14.92
                       LMS
                             0.1000
                                           1.0
           17.72
                                           1.0
      3
                             0.0500
                       LMS
      4
           23.04
                       LMS
                             0.0125
                                           1.0
      5
           30.00
                                          30.0
                    BinMel
                                NaN
```

33.96	LMS	0.0125	1.0
49.22	LMS	0.0500	1.0
51.92	LMS	0.1000	1.0
55.72	LMS	0.0250	1.0
60.00	LowMel	NaN	30.0
65.08	LMS	0.0250	1.0
71.18	LMS	0.0500	1.0
75.82	LMS	0.0125	1.0
82.98	LMS	0.1000	1.0
90.00	${\tt MonLMel}$	NaN	30.0
92.54	LMS	0.0500	1.0
100.50	LMS	0.1000	1.0
106.68	LMS	0.0125	1.0
109.10	LMS	0.0250	1.0
120.00	LowMel	NaN	30.0
123.84	LMS	0.0500	1.0
133.44	LMS	0.0125	1.0
137.34	LMS	0.0250	1.0
144.20	LMS	0.1000	1.0
150.00	${\tt MonRMel}$	NaN	30.0
153.68	LMS	0.0250	1.0
156.84	LMS	0.0500	1.0
160.16	LMS	0.1000	1.0
164.32	LMS	0.0125	1.0
180.00	LowMel	NaN	30.0
187.34	LMS	0.0500	1.0
193.32	LMS	0.1000	1.0
198.38	LMS	0.0125	1.0
205.22	LMS	0.0250	1.0
210.00	MonLMel	NaN	30.0
217.28	LMS	0.0125	1.0
219.46	LMS	0.0500	1.0
222.04	LMS	0.1000	1.0
235.78	LMS	0.0250	1.0
240.00	LowMel	NaN	30.0
247.58	LMS	0.1000	1.0
252.50	LMS	0.0125	1.0
256.40	LMS	0.0500	1.0
	LMS	0.0250	1.0
270.00	BinMel	NaN	30.0
275.34	LMS	0.0500	1.0
283.12	LMS	0.0125	1.0
290.70	LMS	0.1000	1.0
295.34	LMS	0.0250	1.0
300.00	LowMel	NaN	30.0
303.06	LMS	0.0500	1.0
309.02	LMS	0.0250	1.0
	49.22 51.92 55.72 60.00 65.08 71.18 75.82 82.98 90.00 92.54 100.50 106.68 109.10 120.00 123.84 133.44 137.34 144.20 150.00 153.68 156.84 160.16 164.32 180.00 187.34 193.32 198.38 205.22 210.00 217.28 219.46 222.04 235.78 240.00 247.58 252.50 256.40 259.24 270.00 275.34 283.12 290.70 295.34 300.00 303.06	49.22       LMS         51.92       LMS         55.72       LMS         60.00       LowMe1         65.08       LMS         71.18       LMS         75.82       LMS         82.98       LMS         90.00       MonLMe1         92.54       LMS         100.50       LMS         106.68       LMS         109.10       LMS         120.00       LowMe1         123.84       LMS         137.34       LMS         137.34       LMS         150.00       MonRMe1         153.68       LMS         156.84       LMS         160.16       LMS         164.32       LMS         193.32       LMS         198.38       LMS         205.22       LMS         210.00       MonLMe1         217.28       LMS         222.04       LMS         235.78       LMS         240.00       LowMe1         247.58       LMS         259.24       LMS         259.34       LMS         295.34       LM	49.22         LMS         0.1000           51.92         LMS         0.1000           55.72         LMS         0.0250           60.00         LowMel         NaN           65.08         LMS         0.0500           71.18         LMS         0.0500           75.82         LMS         0.1000           90.00         MonLMel         NaN           92.54         LMS         0.0500           100.50         LMS         0.0125           109.10         LMS         0.0250           120.00         LowMel         NaN           123.84         LMS         0.0500           133.44         LMS         0.0250           144.20         LMS         0.1000           150.00         MonRMel         NaN           153.68         LMS         0.0250           156.84         LMS         0.0500           160.16         LMS         0.0125           180.00         LowMel         NaN           187.34         LMS         0.0500           193.32         LMS         0.0125           205.22         LMS         0.0500           219.46

53	312.22	LMS	0.1000	1.0
54	315.62	LMS	0.0125	1.0
55	330.00	${\tt MonRMel}$	NaN	30.0
56	335.16	LMS	0.0250	1.0
57	344.96	LMS	0.0500	1.0
58	348.34	LMS	0.0125	1.0
59	351.88	LMS	0.1000	1.0

# 2 Finding STLAB settings


The following process seems to work well for realising the above design.

- 1. Start with the dimmer device (STLAB 2, right eye). That way, we should always be able to optimise the settings to the other device.
- 2. Search for a pair of spectra that maximise melanopsin contrast (assuming standard 32-year-old observer with 10-degree field size), but constraining every other primary to roughly half of its maximum input. This should help ensure a rod saturating baseline and reserve room in the gamut for LMS modulations.
- 3. Create a sub-problem with the constrained primaries and use this to implement the LMS modulations algebraically.
- 4. Optimise all spectra to the other device.
- 5. Tune the spectra to individual observers (optional, probably not required).

```
[2]: import os
     import os.path as op
     import sys
     import numpy as np
     import pandas as pd
     import matplotlib.pyplot as plt
     plt.style.use('bmh')
     from pysilsub.problems import SilentSubstitutionProblem
     from pysilsub.binocular import BinocularStimulationDevice
     S2 = SilentSubstitutionProblem.from_package_data('STLAB_2_York')
     # Primary bounds to use in optimisation
     bounds = [(0.0, 1.0) for primary in range(S2.nprimaries)]
     # Constrain every other primary to half-maximum. Later we will use
     # these primaries as the background spectrum for LMS modulations.
     constrain = [1,3,5,7,9]
     for led in constrain:
         bounds[led] = (.49,.51)
     # Define the problems
```

```
S2.ignore = ['rh']
S2.target = ['mel']
S2.silence = ['sc', 'mc', 'lc']
S2.bounds = bounds
S2.target_contrast = 'max'
# Find a High-Mel solution
solution_2 = S2.optim_solve(global_search=True, niter=10)
# The low and high melanopsin spectrum settings
s2 low mel = solution 2.x[0:10]
s2_high_mel = solution_2.x[10:]
# Assign the background
S2.background = s2_low_mel
# Plot the solutions
fig2 = S2.plot_solution(s2_high_mel)
# These are the contrasts in each eye
print('STLAB_2 photoreceptor contrasts:')
print(S2.print_photoreceptor_contrasts(s2_high_mel))
# Settings for device 2
print(f'STLAB_2 High-Mel settings: {S2.w2s(s2_high_mel)}')
print(f'STLAB 2 Low-Mel settings: {S2.w2s(s2 low mel)}')
----- optim solve ------
> No background specified, will optimise background.
> Aiming to maximise contrast.
> Performing global optimization with basinhopping and SLSQP
warning: basinhopping: local minimization failure
basinhopping step 0: f -0.0762448
basinhopping step 1: f -0.0961016 trial f -0.0961016 accepted 1 lowest f
-0.0961016
found new global minimum on step 1 with function value -0.0961016
warning: basinhopping: local minimization failure
basinhopping step 2: f -0.0688646 trial f -0.0688646 accepted 1 lowest f
-0.0961016
warning: basinhopping: local minimization failure
basinhopping step 3: f -0.0130398 trial_f -0.0130398 accepted 1 lowest_f
-0.0961016
basinhopping step 4: f -0.0443477 trial_f -0.0443477 accepted 1 lowest_f
-0.0961016
warning: basinhopping: local minimization failure
basinhopping step 5: f -0.0190076 trial_f -0.0190076 accepted 1 lowest_f
-0.0961016
```

```
basinhopping step 6: f -0.3004 trial f -0.3004 accepted 1 lowest f -0.3004
found new global minimum on step 6 with function value -0.3004
basinhopping step 7: f -0.40124 trial f -0.40124 accepted 1 lowest f -0.40124
found new global minimum on step 7 with function value -0.40124
warning: basinhopping: local minimization failure
basinhopping step 8: f -0.40124 trial_f -0.293522 accepted 0 lowest_f -0.40124
basinhopping step 9: f -0.404753 trial f -0.404753 accepted 1 lowest f
-0.404753
found new global minimum on step 9 with function value -0.404753
warning: basinhopping: local minimization failure
basinhopping step 10: f -0.406498 trial f -0.406498 accepted 1 lowest f
-0.406498
found new global minimum on step 10 with function value -0.406498
STLAB_2 photoreceptor contrasts:
      -0.000002
sc
       0.000014
mс
1c
      -0.000000
       0.280079
rh
       0.406498
mel
dtype: float64
None
STLAB 2 High-Mel settings: [1011, 2080, 4029, 2080, 3609, 2006, 218, 2088, 4029,
2088]
STLAB 2 Low-Mel settings: [4030, 2050, 0, 2006, 335, 2088, 1990, 2006, 0, 2006]
```



### 2.1 Matching the low and high melanopsin spectra

```
[8]: S1 = SilentSubstitutionProblem.from_package_data('STLAB_1_York')
S2 = SilentSubstitutionProblem.from_package_data('STLAB_2_York')

# Match the spectra in the other device
Sbin = BinocularStimulationDevice(S1, S2)
Sbin.anchor = 'right' # The less bright of the two
Sbin.optim = 'left'
```

```
s2_settings = [s2_low_mel, s2_high_mel]
s1_settings = Sbin.optimise_settings(s2_settings)


# Plot the solutions
fig, axs = plt.subplots(1, 2, figsize=(7.08, 3.54))

# Low mel
S1.predict_multiprimary_spd(s1_settings[0]).plot(ax=axs[0], label=S1.name)
S2.predict_multiprimary_spd(s2_settings[0]).plot(ax=axs[0], label=S2.name)

# High mel
S1.predict_multiprimary_spd(s1_settings[1]).plot(ax=axs[1], label=S1.name)
S2.predict_multiprimary_spd(s2_settings[1]).plot(ax=axs[1], label=S2.name)

axs[0].set_title('Low melanopsin')
axs[1].set_title('High melanopsin')
for ax in axs:
    ax.legend(loc='best')
```


> Optimising settings. Matching left to right



#### 2.2 LMS modulations

To make the LMS modulations, we can create a sub-problem using the primaries that were constrained in the original optimisation.

```
[77]: import copy
     S2_lum = copy.copy(S2)
     S2_{lum.nprimaries} = 5
     S2_lum.primaries = constrain
     S2_lum.calibration = S2_lum.calibration.loc[constrain]
     S2 lum.ignore = ['rh']
     S2_lum.target = ['sc', 'mc', 'lc']
     S2_lum.silence = ['mel']
     S2_lum.background = s2_high_mel[constrain]
     S2_lum.bounds = [(0.0, 1.0) for primary in range(S2_lum.nprimaries)]
     S2_lum.print_problem()
    *********************
    ******* Silent Substitution Problem *********
    *****************
    Device: STLAB 2 (binocular, right eye)
    Observer: ColorimetricObserver(age=32, field_size=10)
    Ignoring: ['rh']
    Silencing: ['mel']
    Targeting: ['sc', 'mc', 'lc']
    Target contrast: None
    0.51
                                                                 ]
[5]: s2_h_lms = S2_lum.make_contrast_modulation(
        frequency=4,
        sampling_frequency=50,
        target_contrast=.15,
        duration=1
     _ = S2_lum.plot_contrast_modulation(s2_h_lms)
```



## 2.3 Matching the LMS modulations

```
[6]: def combine_settings(background, modulation):
    settings = []
    constrain= [1,3,5,7,9]
    for mod in modulation:
        new = background.copy()
        new[constrain] = mod
        settings.append(new)
    return settings

s1_h_lms = Sbin.optimise_settings(combine_settings(s2_high_mel, s2_h_lms))
```

> Optimising settings. Matching left to right

```
[7]: S1.background = s1_settings[1] # High mel
S2.background = s2_settings[1] # High mel

fig, axs = plt.subplots(1, 2, figsize=(7.08, 3.54))
```

```
_ = S1.plot_contrast_modulation(s1_h_lms, ax=axs[0])
_ = S2.plot_contrast_modulation(combine_settings(s2_settings[1], s2_h_lms),
_ ax=axs[1])

fig.suptitle('LMS probes', fontsize=14)
plt.tight_layout()
```

# LMS probes

